



At-Home Testing Kit Instructions



CHECK YOUR KIT FOR THE FOLLOWING



EXTRA URINE
SAMPLE



TRANSPORT TUBE



SPECIMEN BAG



CUP



REQUISITION FORM



LAB PAK/RETURN LABEL

COLLECT & MAIL YOUR SAMPLE



SEE HOW TO USE THE **247LABKIT**



1. Fill the cup with your urine sample.
2. Locate the transport tube. DO NOT POUR THE LIQUID CONTENTS OF THE TRANSPORT TUBE OUT. Use the pipette to transfer urine into the transport tube. Fill the transport tube between the two black lines.
3. Ensure the transport tube and urine cup are secure. Place both the transport tube and the urine cup inside the specimen bag.
4. Complete the requisition form. Ensure all information is correct. Place the requisition form inside the specimen bag.
5. Place your specimen bag with your sample & requisition form inside the 247Labkit. Place your 247Labkit inside the UPS Laboratory Pak. Your kit comes complete with a free economy return shipping label; however, economy shipping can sometimes be unpredictable, and it may take up to 7 days for our lab to receive your sample.